




# LABORATORY LINES

## SHAKERS



**Kalstein**®

A different accompaniment, at your service





# Kalstein®

A different accompaniment, at your service

OUR SERVICES

## Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



### Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



### Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



### #Letsgivemore ❤️

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



### Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



### Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



### Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



## Roller Mixer

### Model Series YR05762



#### Product Description

YR05762 digital roller mixer provides a gentle but highly efficient rocking and rolling motions, ideal for mixing blood samples, viscous substances and liquid/solid suspensions. Used in a variety of applications, prevention of blood coagulation and immune precipitation etc.

#### Features:

- Rocking and rolling motions for complete mixing
- Display show the speed and time
- Speed range of 10-99rpm
- Seven roller design with small footprint and easy clean
- Digital setting speed and time

#### Technical Specifications

|                       |                                   |
|-----------------------|-----------------------------------|
| Model                 | YR05762                           |
| Speed                 | Adjustable, 10-99rpm $\pm$ 10%    |
| Motion                | Simultaneous roll & tilt $\pm$ 5% |
| Operating Temperature | 4°C – 40°C                        |
| Electrical            | 80-230V, 50-60Hz, 50W             |



#### #Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,  
2 Rue Jean Lantier • 75001 Paris •  
+33 1 78 95 87 89 / +33 6 80 76 07 10 •  
<https://kalstein.eu>  
KALSTEIN FRANCE, S. A. S





 **Kalstein**<sup>®</sup>  
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE  
ACQUISITION  
OF A KALSTEIN  
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:

 **Fundación  
Jacinto Convit**  
<https://www.jacintoconvit.org/>


 **ONE  
TREE  
PLANTED**

<https://onetreepanted.org>

 **Humatem**  
Solidarité et coopération  
biomédicale  
<https://www.humatem.org/>

  
**MANIAPURE  
FOUNDATION USA**  
<https://www.maniapurefoundation.org/>

**Articles Kalstein**  
(Social Networks and Official Websites)



**Kalstein Blog**  
Not just be informed, be fully informed  
**The 3 articles of the week!**  
Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>



**1**  
Read this ARTICLE NOW  
**How does the climate affect your microscope?**  
microscopes in a laboratory are one of the most sensitive pieces of equipment that require to maintain in a controlled environment (refrigerated space, without fluctuations and vibrations) to ensure their performance.  
**Kalstein Blog**  
**Article of the week:**  
How does the climate affect your microscope?  
Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>



**2**  
Read this ARTICLE NOW  
**What is the importance of the laboratory refrigerator?**  
Refrigerator in a laboratory is one of the most sensitive pieces of equipment that require to maintain in a controlled environment (refrigerated space, without fluctuations and vibrations) to ensure their performance.  
**Kalstein Blog**  
**Article of the week:**  
What is the importance of the laboratory refrigerator?  
Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>

**Kalstein Training Videos**  
(Youtube Channel and Official Websites)



**Kalstein France - Jacinto Convit Foundation / Allianc...**  
23 views • 3 years ago



**Kalstein's Microscopes**  
445 views • 1 week ago



PARIS - FRANCE  
**HEADQUARTER**  
2 Rue Jean Lantier  
75001 Paris - France

**Tlf:** +33 (0) 01 78 95 87 89  
**Cel:** +33 (0) 6 80 76 07 10  
[sales@kalstein.eu](mailto:sales@kalstein.eu)

<https://kalstein.eu/>  
<https://www.kalstein.fr/>

