



LABORATORY LINES

VISCOMETER



Kalstein®

A different accompaniment, at your service





Kalstein®

A different accompaniment, at your service

OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



Automatic Digital Brookfield Viscometer

Model YR05859 // YR05860



Product Description

- ✓ The range of No. 1, 2, 3, 4, 5, 6, 7spindles: 10-13,000,000 cp
- ✓ Test range from 10 to 13,000,000 mpa.s (cps)
- ✓ High accuracy
- ✓ Automatic detection
- ✓ Easy to operate

Attachments:

1. No. 0, 1, 2, 3, 4, 5, 6, 7 Spindless
2. Adapter
3. Temperature Sensor
4. Protection Frame
5. Foundation Plate
6. Up-Down Frame
7. English Manual
8. Hexagonal Spanner
9. Solid Wrench



#Let'sgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
Optimum Business Center 450 Rue Baden Powell,
• 34000 Montpellier, France.
Tél: +33 467158849 / +33 680760710/ +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A. S



Automatic Digital Brookfield Viscometer

Model YR05859 // YR05860



Technical Specifications

Model	YR05859	YR05860
Measurement Range	10-13,000,000 mpa.s (cp)	
Display	LCD	
Spindles Specification	No. 0, 1, 2, 3, 4, 5, 6, 7 Note: The bigger the number of the Spindles, the bigger the measuring value	
Rotating Speed of Spindles	0.3 / 0.6 / 0.5 / 3 / 6 / 12 / 30 / 60 rpm/min. with 1, 2, 3, 4# Spindles 0.5 / 1 / 2 / 2.5 / 4 / 5 / 10 / 20 / 50 / 100 rpm/min. with 5, 6, 7# Spindles	
Interface	Printer Interface PC Interface Temperature Sensor Interface	
Automatic Switch	Free selection of proper rotating speed or Spindles number.	
Stable Cursor	Strip-type cursor fully charged displayed when reading basically stable	
Measurement Tolerance	± 1%	
Working Condition	Temperature 5-35 degree Relative humidity no more than 80%	
Outline Dimension (mm)	370×325×280 mm	
Packing (mm)	400×360×310 mm	
Net Weight (kg)	6.8 kg	
Gross Weight (kg)	9 kg	
Additional Information	0.3,0.6,1.5,3,6,12,30,60 rpm/min	0.5,1,2,2.5,4,5,10,20,50,100 rpm/min
Mpa.S	10-13,000,000 cps	



#Letsgivingmore ❤️

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
Optimum Business Center 450 Rue Baden Powell,
• 34000 Montpellier, France.
Tit: +33 467158849 / +33 680760710/ +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A. S





 **Kalstein**[®]
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE
ACQUISITION
OF A KALSTEIN
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



**Fundación
Jacinto Convit**

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



Humatem

Solidarité et coopération
biomédicale

<https://www.humatem.org/>



**MANIAPURE
FOUNDATION USA**

<https://www.maniapurefoundation.org/>

Articles Kalstein
(Social Networks and Official Websites)



Kalstein Blog
Not just be informed, be fully informed

The 3 articles of the week!

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



PARIS - FRANCE
HEADQUARTER
2 Rue Jean Lantier
75001 Paris - France

Tlf: +33 (0) 01 78 95 87 89
Cel: +33 (0) 6 80 76 07 10
sales@kalstein.eu

<https://kalstein.eu/>
<https://www.kalstein.fr/>



1
Read this ARTICLE NOW

How does the climate affect your microscope?

microscope is an instrument that allows observing and measuring objects that are too small to be seen with the naked eye. This is achieved through a complex system of lenses and mirrors, which form

Kalstein Blog

Article of the week:
How does the climate affect your microscope?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



2
Read this ARTICLE NOW

What is the importance of the laboratory refrigerator?

refrigerator is a device that is one of the most common pieces of laboratory equipment. Its purpose is to maintain, in a controlled environment (refrigerated space), various fluids and substances

Kalstein Blog

Article of the week:
What is the importance of the laboratory refrigerator?


Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

Kalstein Training Videos
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Alliances

23 views • 3 years ago



Kalstein's Microscopes

445 views • 1 week ago

