



LABORATORY LINES

BIOLOGICAL MICROSCOPE



Kalstein®

A different accompaniment, at your service





Kalstein®

A different accompaniment, at your service

OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



Biological Microscope
Model YR0237



 **Kalstein**
Biological Microscope
YR0237



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
Optimum Business Center 450 Rue Baden Powell,
• 34000 Montpellier, France.
Tit: +33 467158849 / +33 680760710/ +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A. S





Biological Microscope

Model YR0237



Description:

- ✓ Excellent image quality with infinitive optical system.
- ✓ Comfortable operating with ergonomic structure.
- ✓ Better illumination with sliding-in centerable condenser, Convenient replacement with phase contrast and dark field condenser.

Feature:

- ✓ Extra Wide Field Eyepiece EW10X/20 with Diopter Adjustment as Optional.
- ✓ Comfortable Observation with Eyepoint and interpupillary adjustable.
- ✓ Sliding-in Centerable Condense Adjustable.
- ✓ Handle for easy carrying.
- ✓ Lamp Replacement by Direct Removing Housing.

Optical System	Infinite Optical System	•
Viewing Head	Seidentopf Binocular Viewing Head, 30° Inclined, Interpupillary 48-75mm	•
	Seidentopf Trinocular Viewing Head, 30° Inclined, Interpupillary 48-75mm	○
Eyepiece	Wide Field Eyepiece WF 10×/18	•
	Extra Wide Field Eyepiece EW10×/20 with Diopter Adjustment	○
Objective	Infinite Semi-plan Achromatic Objectives 4×, 10×, 40×, 100×	•
	Infinite Plan Achromatic Objectives 4×, 10×, 40×, 100×	○
Nosepiece	Quadruple Nosepiece	•
Stage	Double Layers Mechanical Stage 140mm×140mm/ 75mm×50mm	•
Condenser	Sliding-in Centerable condenser NA1.25	•
Focusing System	Coaxial Coarse and Fine Adjustment, Fine Division 0.002mm, Coarse Stroke 37.7mm per Rotation, Fine Stroke 0.2mm per Rotation, Moving Range 20mm	•
Illumination	6V/ 20W Halogen Lamp, Brightness Adjustable	•
	3W LED Illumination	○
Phase Contrast Kit		○
Fluorescent Attachment	YX-2 Epi-fluorescent Attachment	○
Dark Field Attachment		○

Note: •Standard Outfit, ○Optional



#Let'sgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
Optimum Business Center 450 Rue Baden Powell,
• 34000 Montpellier, France.
Tit: +33 467158849 / +33 680760710/ +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A. S



Biological Microscope

Model YR0237



Trinocular Model (YR0237-1)



Binocular Model (YR0237)



#Letsgivingore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
Optimum Business Center 450 Rue Baden Powell,
• 34000 Montpellier, France.
Tél: +33 467158849 / +33 680760710 / +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A. S





 **Kalstein**[®]
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE
ACQUISITION
OF A KALSTEIN
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



Fundación
Jacinto Convit

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



Humatem

Solidarité et coopération
biomédicale


<https://www.humatem.org/>



**MANIAPURE
FOUNDATION USA**

<https://www.maniapurefoundation.org/>

Articles Kalstein
(Social Networks and Official Websites)



Kalstein Blog
Not just be informed, be fully informed
The 3 articles of the week!
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



1
Read this ARTICLE NOW
How does the climate affect your microscope?
A microscope is an instrument that allows observing and measuring objects that are too small to be seen with the naked eye. This is achieved through a complex system of lenses and mirrors, which form...

Kalstein Blog
Article of the week:
How does the climate affect your microscope?
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



2
Read this ARTICLE NOW
What is the importance of the laboratory refrigerator?
A refrigerator is a device that is one of the most common pieces of laboratory equipment. It is used to maintain, in a controlled environment (refrigerated space), various fluids and substances...

Kalstein Blog
Article of the week:
What is the importance of the laboratory refrigerator?
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

Kalstein Training Videos
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Alliances

23 views • 3 years ago



Kalstein's Microscopes

445 views • 1 week ago



MONTPELLIER - FRANCE HEADQUARTERS

Optimum Business Center,
450 Rue Baden Powell,
34000 Montpellier - France.

Fax: +33 (0) 4 67 15 88 49

Tel: +33 (0) 6 80 76 07 10

sales@kalstein.eu

<https://kalstein.eu/>
<https://www.kalstein.fr/>

Official Websites for Latin America.

