




# MEDICAL SECTOR

## ULTRASOUND SCANNER



**Kalstein**®

A different accompaniment, at your service





# Kalstein®

A different accompaniment, at your service

OUR SERVICES

## Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



### Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



### Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



### #Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



### Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



### Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



### Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



**KALSTEIN® Ultrasound Scanner**  
**Color Doppler**  
**Model YR05150**



**#Letsgivemore** ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,  
2 Rue Jean Lantier • 75001 Paris •  
+33 1 78 95 87 89 / +33 6 80 76 07 10 •  
<https://kalstein.eu>  
**KALSTEIN FRANCE, S. A. S**



## KALSTEIN® Ultrasound Scanner Color Doppler Model YR05150



### Product Description:

- ✓ Large 15-inch color LED display with true Doppler function, USB ports and VGA port and 2 probe connectors.
- ✓ Built-in battery, which can be worked contentiously at least 3 hours when power off.
- ✓ Main Application:  
Abdomen/Cardiac/Obstetrics/Gynecology/Urology/Andrology/ Small Parts/Vascular/Pediatrics/Musculoskeletal and so on
- ✓ PICTURE MODES: B, 2B, 4B, M, B/M, B/C, B/D, B/C/D, B/CFM/D, CF+B model simultaneously, PDI Color, Dual Color, Simultaneous Sim2D / 3D / 4D Color Composite, PW, CW Duplex / Triplex, CFM, CDE, PD, Directional PD, CD. Anatomy M, Color mode M
- ✓ Color Doppler based on PC with Windows system, which can be connected with any printer of any brand. The printed area is adjustable, which can be images, reports, or image + report, etc.
- ✓ Built-in DICOM 3.0 protocol
- ✓ Multi-language function: English, Chinese, Spanish and Portuguese, Russian, Arabic, French, Vietnamese, Indonesian, German, Persian, Thai, etc.

### Technical Description:

<b>Model number</b>	YR05150
<b>Warranty</b>	2 years
<b>After Sales Service</b>	Return and replacement
<b>Instrument Classification</b>	Class II
<b>Type</b>	Good quality portable color doppler portable ultrasonic diagnostic devices
<b>Image Mode</b>	Model 2D, 3D, CF + B simultaneously, PDI, DPDI, TDI, TSI, etc.
<b>Control of user rights</b>	YES
<b>Probe Connectors</b>	two
<b>PW duplex and triplex freely switched</b>	YES
<b>Automatic IMT measurement</b>	YES
<b>CW function</b>	YES
<b>Complex model image</b>	YES
<b>Automatic PW Tracking Measurements</b>	YES
<b>software package</b>	General, OB/GYN. Small Parts, Urology, Cardiac



**#Letsgivemore** ♥  
Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,  
2 Rue Jean Lantier • 75001 Paris •  
+33 1 78 95 87 89 / +33 6 80 76 07 10 •  
<https://kalstein.eu>  
**KALSTEIN FRANCE, S. A. S**





 **Kalstein**<sup>®</sup>  
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE  
ACQUISITION  
OF A KALSTEIN  
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



**Fundación  
Jacinto Convit**

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



**Humatem**

Solidarité et coopération  
biomédicale


<https://www.humatem.org/>



**MANIAPURE  
FOUNDATION USA**

<https://www.maniapurefoundation.org/>

**Articles Kalstein**  
(Social Networks and Official Websites)



**Kalstein Blog**  
Not just be informed, be fully informed

**The 3 articles of the week!**

Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>



**1**  
Read this ARTICLE NOW

**How does the climate affect your microscope?**

**Kalstein Blog**

**Article of the week:**  
How does the climate affect your microscope?

Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>



**2**  
Read this ARTICLE NOW

**What is the importance of the laboratory refrigerator?**

**Kalstein Blog**

**Article of the week:**  
What is the importance of the laboratory refrigerator?

Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>

**Kalstein Training Videos**  
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Allianc...

23 views • 3 years ago



Kalstein's Microscopes

445 views • 1 week ago



PARIS - FRANCE  
**HEADQUARTER**

2 Rue Jean Lantier  
75001 Paris - France

**Tlf:** +33 (0) 01 78 95 87 89

**Cel:** +33 (0) 6 80 76 07 10  
[sales@kalstein.eu](mailto:sales@kalstein.eu)

<https://kalstein.eu/>  
<https://www.kalstein.fr/>

