




MEDICAL SECTOR

PATIENT AND FETAL MONITORS



Kalstein®

A different accompaniment, at your service





Kalstein®

A different accompaniment, at your service

OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



KALSTEIN® Patient and Fetal Monitor
Doppler Ultrasound Fetal Monitor
Model YR02181



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S



KALSTEIN® Patient and Fetal Monitor Doppler Ultrasound Fetal Monitor Model YR02181



Product Description:

10.2 inch large-screen color TFT display, digital display and waveform oscilloscope of twins heart rate (FHR) and uterine contraction pressure (TOCO) value; Directly viewing 1:1 scale of display and print out waveform; 7-wafer impulse ultrasound transducer, configured with special HR self correlation reading algorithm to ensure correct and reliable fetal hear rate (FHR) measurement; Auto power failure protection and reliable rechargeable Li-ion battery can realize continual monitoring function.

Features:

- ✓ 10.2" Inch large-screen color TFT display, digital display and waveform oscilloscope of twins heart rate (FHR) and uterine contraction pressure (TOCO) value;
- ✓ Directly viewing 1:1 scale of display and print out waveform;
- ✓ 7-wafer impulse ultrasound transducer, configured with special HR self correlation reading algorithm to ensure correct and reliable fetal hear rate (FHR) measurement;
- ✓ Auto power failure protection and reliable rechargeable Li-ion battery can realize continual monitoring function.



#Letsgivemore ❤️

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S



KALSTEIN® Patient and Fetal Monitor
Doppler Ultrasound Fetal Monitor
Model YR02181



Technical Description:

Model	YR02181
Fetal heart rate	
Method of measurement	Ultrasound doppler
Measuring range	50-240 BPM
Measurement tolerance	Less than ± 1 BPM
Scope of the fetal heart rate curve	50-240 BPM
Fetal heart rate numerical display oscilloscope	50-240 BPM
Ultrasound working frequency	2 MHz, 2.5 MHz
Ultrasound output intensity	<5 mW / cm ²
FOR CO	
Method of measurement	External measurement
Measuring range	0-100 unit
Numerical display range of uterine contraction pressure	0-100 unit
Non-linearity tolerance	less than 10%
Auto-zero function	Yes
General specifications	
Fetal movement record	manual control button label
	automatic identification of fetal movement
Display	10.2 "color TFT-LCD (800RGB" 480)
Event mark	9 available
Recording mode	Thermal matrix printer
Recording paper	112mm x 100mm Folded Thermal Printing Paper
Paper feed speed	1,2,3cm / min 3 options
Audible and visual fetal heart rate alarm	low limit 120 BPM, upper limit 160 BPM
Alarm activation time	<3 s
Printer failure alarm	Empty paper alarm
Output volume	Speaker 2-3W Adjustable volume
Power supply	AC: 100-240 V 50/60 Hz ± 1 Hz 40 VA
DC	14.4V / 2000mAh rechargeable Li-ion, continuous work> 4h
Security level	Class I, type B
Size	L296mm x W312 mm x H92 mm
Net working capacity	Yes



#Let'sgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
 2 Rue Jean Lantier • 75001 Paris •
 +33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S





 **Kalstein**[®]
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE
ACQUISITION
OF A KALSTEIN
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



**Fundación
Jacinto Convit**

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



Humatem

Solidarité et coopération
biomédicale


<https://www.humatem.org/>



**MANIAPURE
FOUNDATION USA**

<https://www.maniapurefoundation.org/>

Articles Kalstein
(Social Networks and Official Websites)



Kalstein Blog
Not just be informed, be fully informed

The 3 articles of the week!

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



1
Read this ARTICLE NOW

How does the climate affect your microscope?

Kalstein Blog

Article of the week:
How does the climate affect your microscope?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



2
Read this ARTICLE NOW

What is the importance of the laboratory refrigerator?

Kalstein Blog

Article of the week:
What is the importance of the laboratory refrigerator?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

Kalstein Training Videos
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Allianc...

23 views • 3 years ago



Kalstein's Microscopes

445 views • 1 week ago



PARIS - FRANCE
HEADQUARTER
2 Rue Jean Lantier
75001 Paris - France

Tlf: +33 (0) 01 78 95 87 89
Cel: +33 (0) 6 80 76 07 10
sales@kalstein.eu

<https://kalstein.eu/>
<https://www.kalstein.fr/>

