




# LABORATORY SECTOR

## INCUBATORS



**Kalstein**®

A different accompaniment, at your service





# Kalstein®

A different accompaniment, at your service

OUR SERVICES

## Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



### Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



### Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



### #Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



### Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



### Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



### Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



**KALSTEIN® Incubators**  
**Dry Bath Incubator**  
**Model YR05897 // YR05900**



**#Letsgivemore** ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,  
2 Rue Jean Lantier • 75001 Paris •  
+33 1 78 95 87 89 / +33 6 80 76 07 10 •  
<https://kalstein.eu>  
**KALSTEIN FRANCE, S. A. S**



## KALSTEIN® Incubators Dry Bath Incubator Model YR05897 // YR05900




### Product Description:

YR05894 is a dry block heater controlled by microcomputer. It can be widely used in sample preservation and reaction, DNA amplification and pre-denaturation of electrophoresis, serum coagulation, etc. in fields of medicine, chemical industry, food safety, quality control and environment.

### Features:

- ✓ LED display with timer
- ✓ High precision temperature
- ✓ Built-in over-temperature protect
- ✓ Small size with transparent lid
- ✓ Various blocks can protect samples from contamination



#Letsgivemore ❤️

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,  
2 Rue Jean Lantier • 75001 Paris •  
+33 1 78 95 87 89 / +33 6 80 76 07 10 •  
<https://kalstein.eu>  
KALSTEIN FRANCE, S. A. S



## KALSTEIN® Incubators Dry Bath Incubator Model YR05897 // YR05900



### Technical Specifications:

Model	YR05897	YR05898	YR05899	YR05900
<b>Temperature Control Range</b>	RT+5°C ~150°C			
<b>Heating Up Rate</b>	≤ 10M (20°C to 100°C ) ≤ 30M (20°C to 150°C )			
<b>Temperature Stability @40°C</b>	±0.3°C			
<b>Temperature Stability @100°C</b>	±0.5°C			
<b>Temperature Stability @120°C</b>	±1°C			
<b>Display accuracy</b>	0.1°C			
<b>Timer</b>	0~99h59m			
<b>Max.Temperature</b>	150°C			
<b>Power</b>	200W	400W	600W	400W

### Optional Accessories:

Type	Tube diameter	Capacity	size (W×D×H)
YR05900-op1	6mm	48	96.5×76.5×50 mm
YR05900-op2	7mm	48	96.5×76.5×50 mm
YR05900-op3	10mm	24	96.5×76.5×50 mm
YR05900-op4	12mm	24	96.5×76.5×50 mm
YR05900-op5	13mm	24	96.5×76.5×50 mm
YR05900-op6	15mm	12	96.5×76.5×50 mm
YR05900-op7	16mm	12	96.5×76.5×50 mm
YR05900-op8	19mm	12	96.5×76.5×50 mm
YR05900-op9	20mm	12	96.5×76.5×50 mm
YR05900-op10	26mm	8	96.5×76.5×50 mm
YR05900-op11	28mm	6	96.5×76.5×50 mm
YR05900-op12	38mm	3	96.5×76.5×50 mm
YR05900-op13	0.2ml tubes	48	96.5×76.5×50 mm
YR05900-op14	0.5ml tubes	48	96.5×76.5×50 mm
YR05900-op15	1.5ml tubes	24	96.5×76.5×50 mm
YR05900-op16	2ml tubes	24	96.5×76.5×50 mm
YR05897-op1	0.2ml	96	96.5×153×33.5 mm
YR05897-op2	Flat	96 Elisa	96.5×153×22.5 mm



#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,  
2 Rue Jean Lantier • 75001 Paris •  
+33 1 78 95 87 89 / +33 6 80 76 07 10 •  
<https://kalstein.eu>  
KALSTEIN FRANCE, S. A. S





 **Kalstein**<sup>®</sup>  
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE  
ACQUISITION  
OF A KALSTEIN  
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



Fundación  
Jacinto Convit

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



**Humatem**

Solidarité et coopération  
biomédicale


<https://www.humatem.org/>



**MANIAPURE  
FOUNDATION USA**

<https://www.maniapurefoundation.org/>

**Articles Kalstein**  
(Social Networks and Official Websites)



**Kalstein Blog**  
Not just be informed, be fully informed  
**The 3 articles of the week!**  
Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>



**1**  
Read this ARTICLE NOW  
**How does the climate affect your microscope?**  
microscopes in a laboratory are one of the most sensitive pieces of equipment that require to maintain in a controlled environment (refrigerated space), without fluctuations in humidity and temperature.  
**Kalstein Blog**  
**Article of the week:**  
How does the climate affect your microscope?  
Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>



**2**  
Read this ARTICLE NOW  
**What is the importance of the laboratory refrigerator?**  
Refrigerator in a laboratory is one of the most sensitive pieces of equipment that require to maintain in a controlled environment (refrigerated space), without fluctuations in humidity and temperature.  
**Kalstein Blog**  
**Article of the week:**  
What is the importance of the laboratory refrigerator?  
Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>

**Kalstein Training Videos**  
(Youtube Channel and Official Websites)



**Kalstein France - Jacinto Convit Foundation / Alliances**  
23 views • 3 years ago



**Kalstein's Microscopes**  
445 views • 1 week ago



PARIS - FRANCE  
**HEADQUARTER**  
2 Rue Jean Lantier  
75001 Paris - France

**Tlf:** +33 (0) 01 78 95 87 89  
**Cel:** +33 (0) 6 80 76 07 10  
[sales@kalstein.eu](mailto:sales@kalstein.eu)

<https://kalstein.eu/>  
<https://www.kalstein.fr/>

