




LABORATORY LINES

ICE MAKER



Kalstein®

A different accompaniment, at your service





Kalstein®

A different accompaniment, at your service

OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



#Letsgivemore ♥

Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



Flake Ice Maker Machine

Model YR05143



Features:

- ✓ For fruit and fresh foods, Ice Flakes makes all the difference. It looks good, thanks to its versatility, and tastes good, thanks to its freshness.
- ✓ As a practical cooling medium, Ice Flakes also has great potential for use in the pharmaceutical industry and scientific labs in general.
- ✓ The flakes are produced by machines that can be installed anywhere. The mechanism is simple. Costs are low. The machine is also designed for speedy maintenance operations.
- ✓ High grade stainless steel is used which are anti-corrosive and durable. The compact structure is simple and space-saving.
- ✓ The cabinet are filled with fluoride-free foam thermal cover featuring good insulation effect. The inner liner is food grade plastic. They are both energy-saving and environment-friendly.
- ✓ World famous brand CFC free R134a compressors are used. Major parts and components have passed relevant safety certifications. The product has stable quality and long service life.
- ✓ The ice-making process is fully computer controlled, using imported computer chips to make control reliable and operation easily.
- ✓ The advanced technologies and imported accessories makes noise low, operation smooth and reliable. The ventilation holes and fans are deployed to ensure that motor of reducer can work reliably with high ambient temperature.
- ✓ The patented technology of line of chip ice-making chamber separated from the evaporator which provides the high cooling efficiency and large ice-making capacity.
- ✓ The screw extrusion hob-style ice-making type makes a compact structure, and automatic separation of ice from water. The optimal design of skate blade makes the shape of ice small and practical.
- ✓ The unique float-type water tank water inflow system ensures that there is no residual water. Ice clearance process not needed, and no water wasted, saving water and energy.
- ✓ Fully automatic running with ice full display, water shortage display, in operation display and malfunction warning display, etc. for shutdown protection. They will stop automatically in case of ice full or water shortage and start automatically when ice is not full or water supply is resumed.
- ✓ Flake ice cools more quickly than other ice forms, which can infiltrate into narrow space, featuring fast cooling speed, good ice bathing effect.
- ✓ Long-term performance tests and safety electrical tests will be performed before delivery, ensuring excellent performance of the product.



#Letsgivingore ♥

Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S



Flake Ice Maker Machine

Model YR05143



Model:	YR05143	YR05144
Ice Making Capacity: Kg/24h	300Kg/24h	500Kg/24h
Ice Storage:	200Kg	300Kg
Voltage:	220V/50Hz or 110V/60HZ (should pay extra charge)	
Power:	1360W	2100W
Cooling type:	Air Cooling	
Refrigerant:	R404a	
Water Supply:	Tap Water	
Net Weight: Kg	95	125
Gross Weight: Kg (Carton)	105	135
Product Dimension: (mm)	560*800*1700 mm	760*810*1870 mm
Packaging Dimension: (mm)	650*660*760 mm (head)	845*705*840 mm (head)
	650*650*850 mm (ice bin)	800*850*1060 mm (ice bin)
	20'GP	19
Stuffing Qty	40'GP	38
	40'HQ	42



#Letsgivingore ♥

Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S





 **Kalstein**[®]
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE
ACQUISITION
OF A KALSTEIN
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



**Fundación
Jacinto Convit**

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



Humatem

Solidarité et coopération
biomédicale

<https://www.humatem.org/>



**MANIAPURE
FOUNDATION USA**

<https://www.maniapurefoundation.org/>

Articles Kalstein
(Social Networks and Official Websites)



Kalstein Blog
Not just be informed, be fully informed
The 3 articles of the week!
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



1
Read this ARTICLE NOW
How does the climate affect your microscope?
A microscope is an instrument that allows observing and measuring objects at a very small scale. This is achieved through a complex system of lenses and mirrors, which form...

Kalstein Blog
Article of the week:
How does the climate affect your microscope?
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



2
Read this ARTICLE NOW
What is the importance of the laboratory refrigerator?
A refrigerator is a device that is one of the most common pieces of equipment used in laboratories. It is a controlled environment (refrigerated space), various fluids and substances...

Kalstein Blog
Article of the week:
What is the importance of the laboratory refrigerator?
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

Kalstein Training Videos
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Alliances
23 views • 3 years ago



Kalstein's Microscopes
445 views • 1 week ago



PARIS - FRANCE
HEADQUARTER
2 Rue Jean Lantier
75001 Paris - France

Tlf: +33 (0) 01 78 95 87 89
Cel: +33 (0) 6 80 76 07 10
sales@kalstein.eu

<https://kalstein.eu/>
<https://www.kalstein.fr/>

