



LABORATORY SECTOR

INCUBATORS



Kalstein®

A different accompaniment, at your service





Kalstein®

A different accompaniment, at your service

OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



KALSTEIN® Incubators
Heated Lid Dry Bath Incubator
Model YR05889 // YR05890



#Letsgivingore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S



KALSTEIN® Incubators
Heated Lid Dry Bath Incubator
Model YR05889 // YR05890



Product Description:

YR05889 – YR05890 are a Dry Bath Incubator with heated lid controlled by a microcomputer, excellent temperature control due to advanced semiconductor technology (Pettier-Based), which can work as a thermal cycler. DH-100 can heat, DC-100 can heat and cool, they can be used for sample preservation and reaction, DNA amplification and protein denaturation treatment, and etc.

Features:

- ✓ Impact design, small size
- ✓ Heated lid function can be on/off according to needs
- ✓ High precise and good uniformity in temperature, fast heating and cooling rate
- ✓ 5 inch touch screen display temperature and time, easy to edit, heating lid temperature can be individually set/displayed



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
 2 Rue Jean Lantier • 75001 Paris •
 +33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S



KALSTEIN® Incubators Heated Lid Dry Bath Incubator Model YR05889 // YR05890



Technical Description:

Model	YR05889	YR05890
Voltage	AC100~240V, 50/60Hz	
Block Temperature Setting Range	R.T. -25°C ~100°C	R.T. +5°C ~100°C
Lid Temperature Setting Range	R.T. +5°C ~110°C	R.T. +5°C ~105°C
Timer	1s~99min59sec (0 means infinite)	
Block Temperature Accuracy	±0.5°C	
Lid Temperature Accuracy	±1.0°C	
Block Temperature Uniformity	±0.5°C	
Display Accuracy	0.1°C	
Heating Up Rate (R.T. 25°C)	37°C -100°C ≥ 7°C /min	
Cooling Rate	100°C -37°C ≤ 8°C /min	—
Max.Steps	10	
Max.Cycles	99	
Max.Power	200W	250W

Optional Accesories:

Type	Capacity
YR05890-op1	96×0.2ml
YR05890-op2	54×0.5ml
YR05890-op3	35×1.5ml
YR05890-op4	35×2.0ml



#Let'sgivemore ❤️

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S





 **Kalstein**[®]
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE
ACQUISITION
OF A KALSTEIN
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:

 **Fundación
Jacinto Convit**
<https://www.jacintoconvit.org/>




<https://onetreepanted.org>

 **Humatem**
Solidarité et coopération
biomédicale
<https://www.humatem.org/>


**MANIAPURE
FOUNDATION USA**
<https://www.maniapurefoundation.org/>

Articles Kalstein
(Social Networks and Official Websites)



Kalstein Blog
Not just be informed, be fully informed
The 3 articles of the week!
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



1
Read this ARTICLE NOW
How does the climate affect your microscope?
microscopes in a laboratory are one of the most sensitive pieces of equipment that require to maintain in a controlled environment (refrigerated space), without fluctuations in humidity and temperature.
Kalstein Blog
Article of the week:
How does the climate affect your microscope?
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



2
Read this ARTICLE NOW
What is the importance of the laboratory refrigerator?
Refrigerator in a laboratory is one of the most sensitive pieces of equipment that require to maintain in a controlled environment (refrigerated space), without fluctuations in humidity and temperature.
Kalstein Blog
Article of the week:
What is the importance of the laboratory refrigerator?
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

Kalstein Training Videos
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Alliances

23 views • 3 years ago



Kalstein's Microscopes

445 views • 1 week ago



PARIS - FRANCE
HEADQUARTER

2 Rue Jean Lantier
75001 Paris - France

Tlf: +33 (0) 01 78 95 87 89

Cel: +33 (0) 6 80 76 07 10

sales@kalstein.eu

<https://kalstein.eu/>
<https://www.kalstein.fr/>

