




LABORATORY SECTOR

HOTPLATES AND HOTPLATE STIRRERS



Kalstein®

A different accompaniment, at your service





Kalstein®

A different accompaniment, at your service

OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



KALSTEIN® Hotplates and Hotplate Stirrers
Heater and Hotplate Thermostatic Digital Magnetic Stirrer
Model YR02938



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S



KALSTEIN® Hotplates and Hotplate Stirrers
Heater and Hotplate Thermostatic Digital Magnetic Stirrer
Model YR02938



Product Description:

YR Series are new Digital Hotplates laboratory magnetic stirrer mixer with Heater made of Ceramic top plate , PID Temperature control and digital display which offers excellent chemical resistance .

Product Features:

White ceramic coated top plate provide excellent chemical resistance, thermal-proof , and corrosion resistance, bright white color ceramics coated plate enhance sample visibility.

- ✓ Fixed safety circuit of 400 °C
- ✓ Smooth acceleration of stirring speed ensures no spillage.
- ✓ Digital display for temperature reduce, add and set
- ✓ Heat only, stir only or heat and stir together function
- ✓ Electronic speed control for constant speed even during change in load
- ✓ High magnetic adhesion
- ✓ Various Sizes of the top flexible choice. four sizes , black and white ceramic are all available



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
 2 Rue Jean Lantier • 75001 Paris •
 +33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S



KALSTEIN® Hotplates and Hotplate Stirrers

Heater and Hotplate Thermostatic Digital Magnetic Stirrer

Model YR02938



Technical Description:

Model	YR02938
Top Plate Dimension (mm)	180 x 180
Top Plate Dimension(Inch)	7.09×7.09
Temperature Control	PID Auto control .
Top Plate Shape	Round or Square
Top Plate Material	Ceramic
Temperature Range for machine	Ambient to 400 degree.
Max Temperature for top plate	325 degree
Top Material	White ceramics and Black Ceramic are Both available .
Temperature Display (°C)	Digital Display
Temperature control	PID Auto control
Temperature accuracy	
Thermocouple probe /sensor type	PT100 , with inner and outer thermocouple .
Stir Speed	100- 1800 RPM
Stir Display	Digital display
Stir speed accuracy	
Permissible Ambient Temperature	Ambient Temperature to 40°C
Voltage	160 – 240V , or 110- 130V.
Heating Power supply (w)	550
Motor Power supply	40 W.DC 20
Machine power supply	600
Stirring Flask Capacity (ml)	50ML-5L
Overall Dimensions (mm)	
	391 x 245 x 120
(L × W × H)	
Packing Dimension (mm)	440 x 290 x 170
Shipping Weight (kg)	5.1
Shell material	Time molding waterproof and dustproof shell with Anti-static paint



#Letsgivingore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S





 **Kalstein**[®]
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE
ACQUISITION
OF A KALSTEIN
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:

 **Fundación
Jacinto Convit**
<https://www.jacintoconvit.org/>


 **ONE
TREE
PLANTED**

<https://onetreepanted.org>

 **Humatem**
Solidarité et coopération
biomédicale
<https://www.humatem.org/>

 **MANIAPURE
FOUNDATION USA**
<https://www.maniapurefoundation.org/>

Articles Kalstein
(Social Networks and Official Websites)



Kalstein Blog
Not just be informed, be fully informed

The 3 articles of the week!

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

1
Read this ARTICLE NOW



How does the climate affect your microscope?

microscope is an instrument that allows observing and measuring small objects. This is a very sensitive device that requires a controlled environment to maintain its performance.

Kalstein Blog

Article of the week:
How does the climate affect your microscope?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

2
Read this ARTICLE NOW



What is the importance of the laboratory refrigerator?

Refrigerator is a piece of equipment that is one of the most important in a laboratory. It is used to maintain samples in a controlled environment (refrigerated space), without freezing and without drying.

Kalstein Blog

Article of the week:
What is the importance of the laboratory refrigerator?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

Kalstein Training Videos
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Allianc...

23 views • 3 years ago



Kalstein's Microscopes

445 views • 1 week ago



PARIS - FRANCE
HEADQUARTER

2 Rue Jean Lantier
75001 Paris - France

Tlf: +33 (0) 01 78 95 87 89
Cel: +33 (0) 6 80 76 07 10
sales@kalstein.eu

<https://kalstein.eu/>
<https://www.kalstein.fr/>

