




LABORATORY LINES

HEATING MANTLE



Kalstein®

A different accompaniment, at your service





Kalstein®

A different accompaniment, at your service

OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



KALSTEIN® Heating Mantle**Model YR05052/ YR05053****Product Description**

YR multi heating mantles are designed for repetitive extraction (Kjeldahl, Soxhlet and other extraction procedures), laboratory reflux and distillation applications.

Features

- ✓ Space saving design of the extraction mantle
- ✓ Independent and separate temperature control at each different condition
- ✓ Inner heating element consists soft glass fibre, protected from damage of flask
- ✓ Precise PID temperature control meets high accuracy experiment needs
- ✓ Multipoint heating mantle maximum temperature for heating is 380°C
- ✓ Multipoint heating mantle is widely used for heating in the school laboratories, industrial and mining, food processing, biochemistry, agriculture, petrochemicals engineering and pharmaceutical industry.
- ✓ Shaping shell is made of cold rolled steel sheet with surface spraying treatment.
- ✓ Nichrome heating wire, litter warming up time, long service life.
- ✓ Easy operation, wide application.
- ✓ Precise PID Temperature control meets high accuracy experiments needs.
- ✓ Outer sensor optional.



Model	YR05049	YR05050	YR05051	YR05052	YR05053	YR05054
Capacity (ml)	3 x 250 ml	3 x 500 ml	3 x 1 000ml	6 x 250 ml	6 x 500 ml	6 x 1 000ml
Voltage (V)	200 – 240V / 50- 60Hz or 100 – 120V / 50 – 60Hz					
Max use Temperature (°C)	380					
Power (W)	150x3	250x3	350x3	150x6	250x6	350x6
Working Time	Continuous					
Exterior Size	500 x 290 x 120	650x350x120	900x290x120	1 200x350x120		
Packing Size (mm)	550 x 440 x 170	700x500x170	950x440x170	1 250x500x170		
Net Weight (Kg)	8	14	8	14		
Gross Weight (Kg)	10	16	10	16		

**#Letsgivemore**

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S





 **Kalstein**[®]
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE
ACQUISITION
OF A KALSTEIN
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



**Fundación
Jacinto Convit**

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



Humatem

Solidarité et coopération
biomédicale


<https://www.humatem.org/>



**MANIAPURE
FOUNDATION USA**

<https://www.maniapurefoundation.org/>

Articles Kalstein
(Social Networks and Official Websites)



Kalstein Blog
Not just be informed, be fully informed
The 3 articles of the week!
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



1
Read this ARTICLE NOW
How does the climate affect your microscope?
A microscope is an instrument that allows observing and measuring objects at a very small scale. This is achieved through a complex system of lenses and mirrors, which form...

Kalstein Blog
Article of the week:
How does the climate affect your microscope?
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>




2
Read this ARTICLE NOW
What is the importance of the laboratory refrigerator?
A refrigerator is a device that is one of the most common pieces of equipment used in laboratories. It is a controlled environment (refrigerated space), various fluids and substances...

Kalstein Blog
Article of the week:
What is the importance of the laboratory refrigerator?
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

Kalstein Training Videos
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Alliances
23 views • 3 years ago



Kalstein's Microscopes
445 views • 1 week ago



PARIS - FRANCE
HEADQUARTER
2 Rue Jean Lantier
75001 Paris - France

Tlf: +33 (0) 01 78 95 87 89
Cel: +33 (0) 6 80 76 07 10
sales@kalstein.eu

<https://kalstein.eu/>
<https://www.kalstein.fr/>

