



LABORATORY LINES

CENTRIFUGES



Kalstein®

A different accompaniment, at your service





Kalstein®

A different accompaniment, at your service

OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



KALSTEIN® High-speed table centrifuge

Model YR137-2



Application





















Kalstein Table-top high speed centrifuge has a large capacity. It is widely used in the field of clinical medicine, biochemistry, genetic engineering and immunology, etc. It is essential equipment used by hospitals, scientific research units and institutions of higher learning.

Product Description

- ✓ This multipurpose centrifuge has advanced and intelligent technology.
- ✓ Control by microcomputer, touch screen LCD.
- ✓ It has 10 acceleration and deceleration level files for your choice.
- ✓ The new electronic lock guarantees the security rotation. It has speed speeding device.
- ✓ The new electronic lock guarantees the security rotation.
- ✓ There is an excess speed safety device, an excess temperature safety device and an automatic alarm device.
- ✓ Silastic sealed ring comfort according to GMP. Certified by the US FDA UU



Accessories

<p>No.4</p>  <p>Rotor capacity: 500mlx4 Max Speed: 4000r/min Max RCF: 3220xg</p> 		<p>Name: Microplate bucket Capacity: 4X2X96 well</p>		
	 <p>Bucket</p>	 <p>1.5mlx25x4</p>	 <p>5mlx28x4 Vacuum Tube</p>	 <p>50mlx4x4</p>
	 <p>500ml Bucket</p>	 <p>10mlx20x4</p>	 <p>15mlx10x4</p>	 <p>250mlx4</p>
	 <p>500ml Bottle</p>	 <p>5ml/7mlx20x4 Vacuum Tube</p>	 <p>15mlx8x4</p>	 <p>20mlx9x4</p>
	 <p>4x24x5ml Vacuum Tube</p>	 <p>50mlx4x4</p>	 <p>100mlx2x4</p>	 <p>300mlx4</p>
	 <p>4x37x7ml GRIA Tube</p>			



#Letsgivemore 

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.


















































All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S.



KALSTEIN® High-speed table centrifuge

Model YR137-2



<p>No.1</p>  <p>Rotor Capacity: 1.5ml/2.2mlx12 Max Speed: 20500r/min Max RCF: 29200xg Adapters: 0.2ml, 0.5ml</p>	<p>No.2</p>  <p>Rotor Capacity: 10mlx12 Max Speed: 15000r/min Max RCF: 23120xg Adapters: 1.5ml, 5ml</p>	<p>No.3</p>  <p>Rotor Capacity: 50mlx6 Max Speed: 13000r/min Max RCF: 17940xg Adapters: 15ml, 20ml</p>	<p>No.7</p>  <p>Rotor Capacity: 1.5ml/2.2mlx24 Max Speed: 16500r/min Max RCF: 24890xg Adapters: 0.2ml, 0.5ml</p>	<p>No.8</p>  <p>Rotor Capacity: 100mlx4 Max Speed: 11000r/min Max RCF: 11720xg Adapters: 10ml, 20ml</p>
<p>No.4</p>  <p>Rotor capacity: 750mlx4 Max Speed: 4000r/min Max RCF: 3780xg</p> 	<p>750ml Bucket</p>  <p>750ml Bottle</p> 	<p>5ml / 7mlx24x4 Blood Vacuum Tube</p>  <p>15mlx14x4</p>  <p>50mlx5x4</p>  <p>300mlx4</p> 	<p>5ml / 7mlx37x4 GRIA Tube</p>  <p>20mlx12x4</p>  <p>100mlx3x4</p>  <p>500mlx4</p> 	
<p>No.5</p>  <p>Rotor capacity: 250mlx6 Max Speed: 4200r/min Max RCF: 3500xg</p> 	<p>250mlx6</p>  <p>250ml Bottle</p> 	<p>1.5mlx9x6</p>  <p>5/7mlx9x6 Blood Vacuum Tube</p>  <p>15mlx5x6</p>  <p>50mlx6</p> 	<p>5mlx7x6</p>  <p>10mlx7x6 Internal thread</p>  <p>20mlx4x6</p>  <p>100mlx6</p> 	
<p>No.6</p>  <p>Rotor capacity: 250mlx4 Max Speed: 4300r/min Max RCF: 3582xg</p>	<p>250mlx4</p>  <p>250ml Bottle</p> 	<p>1.5mlx9x4</p>  <p>5/7mlx9x4 Blood Vacuum Tube</p>  <p>15mlx5x4</p>  <p>50mlx4</p> 	<p>5mlx7x4</p>  <p>10mlx7x4 Internal thread</p>  <p>20mlx4x4</p>  <p>100mlx4</p> 	
<p>No.4</p>  <p>Rotor Capacity: 4x2x96 well Max Speed: 4000r/min Max RCF: 3040xg</p>	<p>4x2x96 well Microplate Bucket</p> 	<p>1.5mlx25x4</p>  <p>15mlx10x4</p> 	<p>5mlx28x4 Vacuum tube</p>  <p>50mlx4x4</p> 	<p>10mlx20x4 Internal Thread</p>  <p>It is customizable</p>



#Letsgivemore ❤️

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S.A. S



KALSTEIN® High-speed table centrifuge

Model YR137-2



Model	YR137-2
Maximum speed	5,800r/min
RCF Max	5 060 x g
Maximum capacity	4 x 750ml
Speed Accuracy	±50r/min
Time Settings	1 min to 99h 59min
Noise	<65dB(A)
Power supply	AC 220V±22V 50/60Hz 10A
Dimension (W×D×H)	545 x 660 x 490 (mm)
Packing Size (L x W x H)	670 x 780 x 630 (mm)
Net weight	72Kg

**#Letsgivemore** ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S





A different accompaniment, at your service

WITH THE ACQUISITION OF A KALSTEIN EQUIPMENT

YOU MAKE A CONTRIBUTION TO:

#LETSGIVEMORE 



Fundación Jacinto Convit

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



Humatem

Solidarité et coopération biomédicale

<https://www.humatem.org/>




MANIAPURE FOUNDATION USA

<https://www.maniapurefoundation.org/>

Thank you!

Articles Kalstein
(Social Networks and Official Websites)



Kalstein Blog
Not just be informed, be fully informed

The 3 articles of the week!

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



1
Read this ARTICLE NOW

How does the climate affect your microscope?

microscope is an instrument that allows observing and measuring objects that are too small to be seen with the naked eye. This is achieved through a complex system of lenses and mirrors, which form

Kalstein Blog

Article of the week:
How does the climate affect your microscope?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



2
Read this ARTICLE NOW

What is the importance of the laboratory refrigerator?

refrigerator is a device that is one of the most common pieces of equipment used in laboratories. It is a controlled environment (refrigerated space), various fluids and substances

Kalstein Blog

Article of the week:
What is the importance of the laboratory refrigerator?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

Kalstein Training Videos
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Alliances
23 views • 3 years ago



Kalstein's Microscopes
445 views • 1 week ago



PARIS - FRANCE
HEADQUARTER
2 Rue Jean Lantier
75001 Paris - France

Tlf: +33 (0) 01 78 95 87 89
Cel: +33 (0) 6 80 76 07 10
sales@kalstein.eu

<https://kalstein.eu/>
<https://www.kalstein.fr/>

