



LABORATORY LINES

CHAMBERS



Kalstein®

A different accompaniment, at your service





Kalstein®

A different accompaniment, at your service

OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



Humidity Chamber Model YR05802 // YR05807



Product Description

Without fluorine refrigerant system, reasonable air circulation system to ensure temperature uniformity; intelligent control to ensure constant and precision of temperature and humidity; large capacity humidifying system; liner with lighting device.

With constant temperature and humidity control function, it is a based equipment for environment protection, medical, drug testing, health and epidemic prevention, livestock, aquaculture and other scientific research and college departments.

Recommended the use of biological engineering and water analysis, Industrial performance test field etc..

Features

1. Features
2. The inner chamber is made of mirror stainless steel.
3. Semicircular corners inner chamber is easy to clean.
4. Shelf can be adjusted.
5. LED microcomputer controller make the temperature and humidity constant.
6. Touch key parameter setting.
7. Timer, over temperature alarm function.
8. Large angle observation window, easy to watch the goods in the inner chamber with lighting.
9. Large-capacity humidifier, easy to use

Product introduction

With constant temperature and humidity control function. It is a based equipment for environmental protection, medical, drug testing, health and epidemic prevention, livestock, aquaculture and other scientific research and college departments. Also, widely application in biological engineering and water analysis field etc.



#Let'sgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
Optimum Business Center 450 Rue Baden Powell,
• 34000 Montpellier, France.
Tit: +33 467158849 / +33 680760710/ +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A. S



Humidity Chamber
Model YR05802 // YR05807



Technical Specifications

Model	YR05802	YR05803	YR05804
Temperatur range (°C)	0~60		
Accuracy (°C)	±0.1		
Timing range	0~999 min		
Fluctuation (%)	±0.5		
Voltage (V)	220		
Power (KW)	0.4	0.6	0.8
Inner size (mm)(H*W*D)	580×450×380/100L	800×500×380/150L	1000×550×450/250L
Overall dimension (mm) (H*W*D)	1100×590×620	1350×640×660	1550×690×660
Packing size (mm) (H*W*D)	1220×670×700	1470×720×740	1670×770×740
Net weight (KG)	93	105	133
Shipping weight (KG)	110	125	155
Model	YR05805	YR05806	YR05807
Temperatur range (°C)	0~60		
Accuracy (°C)	±0.5		
Timing range			
Fluctuation (%)	±1		
Voltage (V)	220		
Power (KW)	0.4	0.6	0.8
Inner size (mm)(H*W*D)	580×450×380/100L	800×500×380/150L	1000×550×450/250L
Overall dimension (mm) (H*W*D)	1100×590×620	1350×640×660	1550×690×660
Packing size (mm) (H*W*D)	1220×670×700	1470×720×740	1670×770×740
Net weight (KG)	100	110	140
Shipping weight (KG)	120	140	165



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
 Optimum Business Center 450 Rue Baden Powell,
 • 34000 Montpellier, France.
 Tit: +33 467158849 / +33 680760710/ +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A. S



Humidity Chamber Model YR05802 // YR05807



Printer



Communication Software



Test hole



#Letsgivemore ♥


Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
Optimum Business Center 450 Rue Baden Powell,
• 34000 Montpellier, France.
Tél. +33 467158849 / +33 680760710 / +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A. S




Humidity Chamber Model YR05802 // YR05807



Built-in large capacity humidification system

The built-in large-capacity humidification system automatically fills the water to ensure stable fluctuations.



Tempered glass door

Built-in toughened glass door, clear observation of window items



Upgraded LCD display

Programmable microcomputer technology with timing over-temperature alarm




Fluorine-free balanced refrigeration system

Fluorine-free refrigeration, energy-efficient, with inlet and outlet



Lighting device

Hood lighting device for easy observation



Air circulation system

Air circulation system ensures temperature and humidity stability in the working room.



#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S., Optimum Business Center 450 Rue Baden Powell, • 34000 Montpellier, France.
Tél: +33 467158849 / +33 680760710 / +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A. S





 **Kalstein**[®]
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE
ACQUISITION
OF A KALSTEIN
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:

 **Fundación
Jacinto Convit**
<https://www.jacintoconvit.org/>


 **ONE
TREE
PLANTED**

<https://onetreepanted.org>

 **Humatem**
Solidarité et coopération
biomédicale
<https://www.humatem.org/>

 **MANIAPURE
FOUNDATION USA**
<https://www.maniapurefoundation.org/>

Articles Kalstein
(Social Networks and Official Websites)



Kalstein Blog
Not just be informed, be fully informed
The 3 articles of the week!
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



1
Read this ARTICLE NOW
How does the climate affect your microscope?
A microscope is an instrument that allows observing and measuring objects that are too small to be seen with the naked eye. This is achieved through a complex system of lenses and mirrors, which form
Kalstein Blog
Article of the week:
How does the climate affect your microscope?
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>




2
Read this ARTICLE NOW
What is the importance of the laboratory refrigerator?
A refrigerator is a device that is one of the most common pieces of equipment used in laboratories to maintain, in a controlled environment (refrigerated space), various fluids and substances.
Kalstein Blog
Article of the week:
What is the importance of the laboratory refrigerator?
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

Kalstein Training Videos
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Alliances
23 views • 3 years ago



Kalstein's Microscopes
445 views • 1 week ago



PARIS - FRANCE
HEADQUARTER
2 Rue Jean Lantier
75001 Paris - France

Tlf: +33 (0) 01 78 95 87 89
Cel: +33 (0) 6 80 76 07 10
sales@kalstein.eu

<https://kalstein.eu/>
<https://www.kalstein.fr/>

