




MEDICAL SECTOR

INFANT PHOTOTHERAPY INCUBATOR



Kalstein®

A different accompaniment, at your service





Kalstein®

A different accompaniment, at your service

OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.




Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



KALSTEIN® Infant Phototherapy Incubator
Infant Phototherapy Incubator
Model YR02182



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S.



KALSTEIN® Infant Phototherapy Incubator Infant Phototherapy Incubator

Model YR02182



Product Description:

The YR02182 mode is the basic one of YR02182 series with the function of infant incubator and phototherapy treatment for neonate patients with jaundice. In addition, the unit can show the accurate parameters, air temperature, skin temperature and multi alarms with three levels.

Application:

The unit combines the incubator with phototherapy treatment for neonate patient that using in NICU. We, Kalstein also have patient monitor for newborn using in NICU too. Widely applicable in all hospitals and clinics and optimize for preterm births or for some ill full-term babies or jaundiced baby.

Features:

- ✓ Phototherapy light intensity is constantly adjustable
- ✓ Time of bilirubin phototherapy can be paused, recorded and accumulated
- ✓ Phototherapy treatment duration can be pre-set and it automatically shuts off when a pre-set duration is reached
- ✓ Working life of LED phototherapy unit is lasting 20,000 hours
- ✓ Silicone mattress protect neonate from came while phototherapy treatment
- ✓ Microprocessor servo-controller ensures a precise and reliable temperature control
- ✓ Servo-controlled humidification system maintains the relative humidity at a constant level in the incubator
- ✓ Set temperature, monitored air temperature and baby skin temperature and heating power are distinguished in display
- ✓ >37°COverride temperature function
- ✓ Self-checking function for safety affirmative
- ✓ Audible and visible alarm indications for multiple failures with 3 priorities
- ✓ Chord alarm notes nurse farther than 8meters, but leaves asleep baby uninterrupted
- ✓ Aluminous framework for constant temperature keeping
- ✓ Continuous adjustment of mattress inclination up to 10°
- ✓ Drawer-type water reservoir is easily removable for cleaning and disinfection
- ✓ Nature air circulation system
- ✓ Cupboard cabinet
- ✓ RS-232 Connector.

Optional Accessories:

- ✓ Cabinet
- ✓ Side door, front-rear door
- ✓ Humidity monitor and controller
- ✓ Double wall hood
- ✓ Vertical height adjustable
- ✓ Infant phototherapy unit (OSRAM tubes)



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S



KALSTEIN® Infant Phototherapy Incubator
Infant Phototherapy Incubator
Model YR02182



Technical Description:

Model	YR02182
Air temperature control	•
Control of skin temperature	•
37 ° air temperature control	•
The baby bed can be tilted	•
Humidity display	•
Humidity control	○
Superior Phototherapy Unit	•
Lower Phototherapy Unit	•
Side door	○
Aluminum tank for heat preservation	•
Windows and access ports operation	1 port, 4 windows
Cabinet	○
Adjustable working height	○
Additional rear door	○
Baby scale bed	○
Silent DC motor	○



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
 2 Rue Jean Lantier • 75001 Paris •
 +33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S





 **Kalstein**[®]
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE
ACQUISITION
OF A KALSTEIN
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



**Fundación
Jacinto Convit**

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



Humatem

Solidarité et coopération
biomédicale


<https://www.humatem.org/>



**MANIAPURE
FOUNDATION USA**

<https://www.maniapurefoundation.org/>

Articles Kalstein
(Social Networks and Official Websites)



Kalstein Blog
Not just be informed, be fully informed
The 3 articles of the week!
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



1
Read this ARTICLE NOW
How does the climate affect your microscope?
Microscopes are instruments used for observing microscopic objects. This is a sensitive device that requires a controlled environment to maintain its performance.
Kalstein Blog
Article of the week:
How does the climate affect your microscope?
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



2
Read this ARTICLE NOW
What is the importance of the laboratory refrigerator?
Refrigeration is a key element of the most advanced equipment that is used in laboratories. It is a controlled environment (refrigerated space), without fluctuations in temperature.
Kalstein Blog
Article of the week:
What is the importance of the laboratory refrigerator?
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

Kalstein Training Videos
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Allianc...

23 views • 3 years ago



Kalstein's Microscopes

445 views • 1 week ago



PARIS - FRANCE
HEADQUARTER

2 Rue Jean Lantier
75001 Paris - France

Tlf: +33 (0) 01 78 95 87 89

Cel: +33 (0) 6 80 76 07 10

sales@kalstein.eu

<https://kalstein.eu/>
<https://www.kalstein.fr/>

