




MEDICAL SECTOR

INFANT RADIANT WARMER



Kalstein®

A different accompaniment, at your service





Kalstein®

A different accompaniment, at your service

OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



KALSTEIN® Infant Radiant Warmer
Infant Radiant Warmer
Model YR02191



#Letsgivingore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S



KALSTEIN® Infant Radiant Warmer Infant Radiant Warmer

Model YR02191



Product Description:

The YR02191 unit is our middle design of infant radiant warmer with a unique thermal environment for premature and full-term neonates. Moreover, the resuscitation accessories can be optional by the unit in emergency conditions.

Application:

The YR02191 model provide appropriate pre-heating, servo-controlled temperature and skin temperature functions reducing the care pressure for caregivers and increasing the care level for newborn patients. Widely applicable in all hospitals and clinics and optimize for preterm births or for some ill full-term babies or jaundiced baby.

Optional Accessories:

- ✓ Oxygen therapy system unit, oxygen hood
- ✓ Suction unit
- ✓ Resuscitation bag
- ✓ Infant head fixating
- ✓ Electric height adjustment
- ✓ Phototherapy unit (LED light or OSRAM tubes)

- ✓ Pre-warm, manual control and baby skin temperature mode
- ✓ Microprocessor based servo controlled temperature system with far-infrared heater
- ✓ Temperature deviation can be observed and amended in front panel for error measurement
- ✓ Self-checking function for safety affirmative
- ✓ Lower phototherapy units
- ✓ The inclination of bassinet can be adjusted.
- ✓ Heater head can be adjusted horizontally between 0-90°
- ✓ Foldable transparent protector
- ✓ Skin temperature sensor fall protect function to avoid over warming
- ✓ Audible and visible alarm indications for multiple failures with 3 priorities
- ✓ Convenient X-ray cassette tray under infant bed
- ✓ Big baby bed is convenient for operation
- ✓ 40 % Constant power warming function
- ✓ Double sides phototherapy unit provide treatment for jaundice
- ✓ The vertical height can be adjusted.
- ✓ Oxygen therapy unit, resuscitation balloon, infant head fixed shelf ,oxygen hood, low pressure suction system.
- ✓ APAGER timer
- ✓ Cylindrical chassis
- ✓ RS-232 connector



#Letsgivemore ❤️

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S



KALSTEIN® Infant Radiant Warmer
Infant Radiant Warmer
Model YR02191



Technical Description:

| | |
|------------------------------------|---|
| Model | YR02191 |
| Temperature control mode | Preheating, manual control and skin temperature |
| Rotating radiant heater | ● |
| Wallbed | ● |
| Fault alarm | ● |
| APGAR timer | ● |
| X-ray cassette | ● |
| Bed dimension | big |
| Superior phototherapy units | ○ |
| Lower phototherapy units | ● |
| Adjustable working height | ○ |
| Low pressure suction unit | ○ |
| Oxygen therapy unit | ○ |
| Oxygen hood | ○ |
| Child headrest unit | ○ |



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
 2 Rue Jean Lantier • 75001 Paris •
 +33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S





 **Kalstein**[®]
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE
ACQUISITION
OF A KALSTEIN
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



**Fundación
Jacinto Convit**

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



Humatem

Solidarité et coopération
biomédicale


<https://www.humatem.org/>



**MANIAPURE
FOUNDATION USA**

<https://www.maniapurefoundation.org/>

Articles Kalstein
(Social Networks and Official Websites)



Kalstein Blog
Not just be informed, be fully informed

The 3 articles of the week!

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



1
Read this ARTICLE NOW

How does the climate affect your microscope?

Kalstein Blog

Article of the week:
How does the climate affect your microscope?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



2
Read this ARTICLE NOW

What is the importance of the laboratory refrigerator?

Kalstein Blog

Article of the week:
What is the importance of the laboratory refrigerator?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

Kalstein Training Videos
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Allianc...

23 views • 3 years ago



Kalstein's Microscopes

445 views • 1 week ago



PARIS - FRANCE
HEADQUARTER

2 Rue Jean Lantier
75001 Paris - France

Tlf: +33 (0) 01 78 95 87 89
Cel: +33 (0) 6 80 76 07 10
sales@kalstein.eu

<https://kalstein.eu/>
<https://www.kalstein.fr/>

