



LABORATORY SECTOR

MICROSCOPE



Kalstein®

A different accompaniment, at your service





Kalstein®

A different accompaniment, at your service

OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



KALSTEIN® Microscope Metallurgical Microscope Model YR0256



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S



KALSTEIN® Microscope Metallurgical Microscope Model YR0256



Product Description:

- ✓ Infinity Optical System, Bright Field/Dark Field, DIC, Polarizing
- ✓ Total Magnification 50x-1000x With Infinity Plan BF/DF Objectives
- ✓ Double Layer Mechanical Stage 186*138mm, Moving Range 74*50mm
- ✓ Maximum Sample Height 30mm/50mm, With Specimen Preparation Plate
- ✓ Reflect, Transmit/Reflect 24V 100W Halogen Light Kohler Illumination



Technical Specifications:

Model	YR0256
Standard Function	BF-DF
Upgrade Function	BF-DF-PL-DIC
Features	Infinity BF+DF DIC
Optical System	Infinity Plan
Head Type	Seidentopf Trinocular
Head Inclination	30°
Head Rotation	Head Rotate 360°
Head Interpupillary	40-75mm
Head Diopter	Left Tube Adjustable
Head Light Path	E50 : P50
Eyepiece Diameter	30.0mm
Nosepiece Features	Inward with DIC Slot
Nosepiece Holes	5 holes
Objective Features	LWD Covreless BF+DF DIC
Objective Type	Infinity Plan Dark Field Metallurgical
Working Stage Features	Size ≥ 140 x 180mm
Working Stage Type	2 Layers Mechanical
Focusing Features	Safety Stop Lock Ring Tension Adjustable
Focusing Type	Coaxial Coarse & Fine
Condenser	Rack & Pinion center adjustable Iris Diaphragm
Features	
Condenser Type	Swing – Out
Condenser N.A	0.9/0.25
Light Source Features	Kohler Illumination Aspheric Collector Halogen 100W
Light Source Type	Halogen
Power Supply Type	AC 110V-240V



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S





 **Kalstein**[®]
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE
ACQUISITION
OF A KALSTEIN
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



**Fundación
Jacinto Convit**

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



Humatem

Solidarité et coopération
biomédicale


<https://www.humatem.org/>



**MANIAPURE
FOUNDATION USA**

<https://www.maniapurefoundation.org/>

Articles Kalstein
(Social Networks and Official Websites)



Kalstein Blog
Not just be informed, be fully informed
The 3 articles of the week!
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



1
Read this ARTICLE NOW
How does the climate affect your microscope?
microscope is an instrument that uses a series of lenses to magnify images of small objects. This is a complex device that requires a lot of maintenance and care to ensure it works properly.

Kalstein Blog
Article of the week:
How does the climate affect your microscope?
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



2
Read this ARTICLE NOW
What is the importance of the laboratory refrigerator?
Refrigerator is a piece of equipment that is used to store and maintain samples in a controlled environment (refrigerated space), without frost and humidity.

Kalstein Blog
Article of the week:
What is the importance of the laboratory refrigerator?
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

Kalstein Training Videos
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Allianc...

23 views • 3 years ago



Kalstein's Microscopes

445 views • 1 week ago



PARIS - FRANCE
HEADQUARTER

2 Rue Jean Lantier
75001 Paris - France

Tlf: +33 (0) 01 78 95 87 89
Cel: +33 (0) 6 80 76 07 10
sales@kalstein.eu

<https://kalstein.eu/>
<https://www.kalstein.fr/>

