


# Kalstein A different accompaniment, at your service

## LABORATORY SECTOR INCUBATORS





All rights reserved ® KALSTEIN France S. A. S., 2 Rue Jean Lantier • 75001 Paris • +33 1 78 95 87 89 /+33 6 80 76 07 10 • https://kalstein.eu KALSTEIN FRANCE, S. A. S





OUR SERVICES

#### **Benefits and Support**

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



### Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



#### **Quick Response**

Our work team is always available to response all your consults or questions, in order to support you in any situation.





#### #Letsgivemore 💚

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



#### **Technical Support**

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.





#### **Delivery Logistics**

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or.air.



#### Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



#### **KALSTEIN®** Incubators

**Mini Dry Bath Incubator** 

Model YR05891 // YR05893







#### #Letsgivemore 💗

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved ® KALSTEIN France S. A. S., 2 Rue Jean Lantier • 75001 Paris • +33 1 78 95 87 89 /+33 6 80 76 07 10 • https://kalstein.eu

KALSTEIN FRANCE, S. A. S



#### **KALSTEIN®** Incubators

**Mini Dry Bath Incubator** 

Model YR05891 // YR05893















#### **Product Description:**

Mini Dry Bath Incubator is controlled by PID, applied to vehicle power supply due to small palm-size and light weight. Widely used in rapid inspection kit or small labs.

#### Features:

- ✓ Display with timer and temperature
- ✓ DC 12V, built-in over temperature protect
- ✓ Automatic failure detection with buzzer alarm
- ✓ Easy to change blocks
- ✓ 9 Stored programs, each contains 3 segments











#Letsgivemore 💗

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved ® KALSTEIN France S. A. S., 2 Rue Jean Lantier • 75001 Paris • +33 1 78 95 87 89 /+33 6 80 76 07 10 • https://kalstein.eu



KALSTEIN FRANCE, S. A. S

#### **KALSTEIN®** Incubators

**Mini Dry Bath Incubator** 

Model YR05891 // YR05893









#### **Technical Specifications**:

Model	YR05891	YR05892	YR05893
Temperature Control	RT+5°C ~100°C	RT-20°C ~100°C	DT - E9C 1009C
Range	K1+5 C ~100 C	R1-20 C ~100 C	RT+5°C ~100°C
Heating Up Rate(20°C to	≤ 15m		
100°C )			
Temperature Control	±0.3°C		
Accuracy			
Temperature Stability	±0.3°C		
@40°C			
Temperature Stability	±0.3°C		
@100°C			
Display Accuracy	0.1°C		
Timer	0~99h, 0~99m, 0~99s		
Max. Temperature	100°C		
Max. Power	48W	60W	120W
Cooling Down Rate	Natural cooling		≤ 30m 100-25°C ≤ 30m RT 20°C
Heating Lid	No		Yes

#### **Optional Accesories:**

Туре	Tube Type	Description	Type Capacity	
YR05893-op1	0.2ml tube	32 x 0.2ml Centrifuge Tube	32	
YR05893-op2	0.5ml tube	24 x 0.5 ml Centrifuge Tube	24	
YR05893-op3	1.5ml tube	15 x 1.5 ml Centrifuge Tube	12	
YR05893-op4	2.0ml tube	15 x2.0 ml Centrifuge Tube	15	
YR05893-op5	Cuette	8~12.5 x 12.5 mm (Cuette)	8	
YR05893-op6	15ml tube (Working without lid)	6 x 15 ml Centrifuge Tube – Open lid when working	6	
YR05893-op7	50ml tube (Working without lid)	2 x 50 ml Centrifuge Tube – Open lid when working	2	
YR05893-op8	5ml tube (Working without lid)	12 x 5.0ml Centrifuge Tube	12	
YR05893-op9		Customized		



#Letsgivemore 💗

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved ® KALSTEIN France S. A. S., 2 Rue Jean Lantier • 75001 Paris • +33 1 78 95 87 89 /+33 6 80 76 07 10 • https://kalstein.eu



KALSTEIN FRANCE, S. A. S





#LETSGIVEMORE >

# WITH THE ADCQUISITION OF A KALSTEIN EQUIPMENT

YOU MAKE A CONTRIBUTION TO:





https://onetreeplanted.org



Solidarité et coopération biomédicale

https://www.humatem.org/



https://www.maniapurefoundation.org/




### ANY DOUBT? Contact us!

#### **Articles Kalstein**

(Social Networks and Official Websites)







Article of the week: How does the climate affect your microscope?

Read more on Kalstein's social medias or official websites
https://kalstein.eu/articles





Article of the week: What is the importance of the laboratory refrigerator?

Read more on Kalstein's social medias or official websites https://kalstein.eu/articles



PARIS - FRANCE HEADQUARTER

2 Rue Jean Lantier 75001 Paris - France

**Tlf:** +33 (0) 01 78 95 87 89 **Cel:** +33 (0) 6 80 76 07 10 sales@kalstein.eu

https://kalstein.eu/ https://www.kalstein.fr/




(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Allianc...

23 views • 3 years ago



Kalstein's Microscopes

445 views • 1 week ago

