



LABORATORY SECTOR

BALANCES



Kalstein®

A different accompaniment, at your service





Kalstein®

A different accompaniment, at your service

OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



KALSTEIN® Balances
Analytical Balance
Model YR05620 // YR05622



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S



KALSTEIN® Balances

Analytical Balance

Model YR05620 // YR05622



Product Description:

- ✓ Appearance design and utility model access to patents
- ✓ 7 inch touch screen
- ✓ Large range and high precision
- ✓ Large size stainless steel pan
- ✓ AC and DC power supply
- ✓ Date/Time set
- ✓ Tare function/Counting
- ✓ Check weighed/Percentage/Density
- ✓ Net and Gross/Sum
- ✓ Unit conversion(g/ct/oz)...
- ✓ Overload alarm
- ✓ Connect PC and printer
- ✓ Optional label printer
- ✓ Function Expansion: Under weighing / Remote display/Internal calibration



Technical Specifications:

Model	YR05620	YR05621	YR05622
Capacity (g)	3000	4000	5000
Sensor	Load Cell		
MIN weighing	0.01		
Stable time	≤3S		
Display	7 inch touch screen		
Pan size	305×230mm		
Operation temp	5-25°C		
Repeat ability	±3d		
Liner	±3d		
Interface	RS232/RJ45/USB...		
In/output	excel, (USB)		
Internal calibration	Internal Calibration (Option)		
Weighing mode	Counting, Checking, Percentage, Density		
Function expansion	Internal Calibration, Under weighing, Remote Display		
Printer	Label printing, label format, content setting, etc.		
Host dimension	395 x 118 x 320mm		
Packing size	415 × 165 × 430mm		
Machine	Input: 220V 50HZ/110V 60HZ		
	Output:12V		
Battery	6V/4Ah Rechargeable Battery		



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S





 **Kalstein**[®]
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE
ACQUISITION
OF A KALSTEIN
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



**Fundación
Jacinto Convit**

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



Humatem

Solidarité et coopération
biomédicale

<https://www.humatem.org/>



**MANIAPURE
FOUNDATION USA**

<https://www.maniapurefoundation.org/>

Articles Kalstein
(Social Networks and Official Websites)



Kalstein Blog
Not just be informed, be fully informed

The 3 articles of the week!

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



1
Read this ARTICLE NOW

How does the climate affect your microscope?

Kalstein Blog

Article of the week:
How does the climate affect your microscope?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



2
Read this ARTICLE NOW

What is the importance of the laboratory refrigerator?

Kalstein Blog

Article of the week:
What is the importance of the laboratory refrigerator?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

Kalstein Training Videos
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Allianc...

23 views • 3 years ago



Kalstein's Microscopes

445 views • 1 week ago



PARIS - FRANCE
HEADQUARTER

2 Rue Jean Lantier
75001 Paris - France

Tlf: +33 (0) 01 78 95 87 89

Cel: +33 (0) 6 80 76 07 10

sales@kalstein.eu

<https://kalstein.eu/>
<https://www.kalstein.fr/>

