

LABORATORY SECTOR

INCUBATORS









our services Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings In any part of the world, receive your

induction or training from our specialized team of engineers



Quick Response Our work team is always available to response all your consults or questions, in order to support you in any situation.





#Letsgivemore 💜

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics We take care of all the necessar

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or.air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.





KALSTEIN® Incubators Shaking Incubator Model YR05759 // YR05759-1

CE





#Letsgivemore 💗

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved ® KALSTEIN France S. A. S., 2 Rue Jean Lantier • 75001 Paris • +33 1 78 95 87 89 / +33 6 80 76 07 10 • https://kalstein.eu KALSTEIN FRANCE, S. A. S



www.kalstein.eu

All rights reserved © KALSTEIN France S. A. S.,

KALSTEIN® Incubators Shaking Incubator Model YR05759 // YR05759-1



Product Description:

Incubator shaker is indispensable laboratory equipment in the plant, biological, microbiological, biological products, genetics, virus, medicine, environmental protection and other scientific research, education and production departments. Universal spring trial bottle rack; Mechanical timer; Steeples speed regulation; Can be equipped with a variety of flasks; test tubes to use; Stable and reliable in continuous work.

Technical Description:

Model	YR05759	YR05759-1
rotation frequency	20-300rpm	
frequency precision	±1rpm	
stroke	24mm	
	500ml*3,	_
standard configuraton	250ml*4,	500ml*28
	100ml*4	
	500ml*9,	1000ml*18, 250ml*43
maximum volume	250ml*16,	
	100ml*25	_
timer range	(0-99h)	
temp control range	RT+5-60°C	
temp resolution	±0.1°C	
temp uniformity	±1°C	
display	LCD	
shaking plate number	1	
Exterior dimensions (mm) (H+W+D)	640*605*520	1110*765*780
Net.Weight (kg)	72	160
power(W)	350	700
power supply	220V±10%	
shaking way	rotating	
motor	brushless motor	



#Letsgivemore 💗

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved ® KALSTEIN France S. A. S., 2 Rue Jean Lantier • 75001 Paris • +33 1 78 95 87 89 / +33 6 80 76 07 10 • https://kalstein.eu KALSTEIN FRANCE, S. A. S







#LETSGIVEMORE

WITH THE **ADCQUISITION OF A KALSTEIN EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:





https://onetreeplanted.org

Humatem Solidarité et coopération biomédicale https://www.humatem.org/



https://www.maniapurefoundation.org/



ANY DOUBT? Contact us!







PARIS - FRANCE **H E A D Q U A R T E R** 2 Rue Jean Lantier 75001 Paris - France

Tlf: +33 (0) 01 78 95 87 89 **Cel:** +33 (0) 6 80 76 07 10 sales@kalstein.eu


https://kalstein.eu/ https://www.kalstein.fr/

Kalstein Training Videos (Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Allianc...

23 views • 3 years ago



Kalstein's Microscopes 445 views • 1 week ago

