




LABORATORY SECTOR

CENTRIFUGES



Kalstein®

A different accompaniment, at your service





Kalstein®

A different accompaniment, at your service

OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



KALSTEIN® Centrifuges
Tabletop High Speed Centrifuge
Model YR019-TG



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S



KALSTEIN® Centrifuges

Tabletop High Speed Centrifuge

Model YR019-TG



Product Description:

- ✓ Microcomputer control, brushless DC engine, stable operation, low noise level and high speed precision.
- ✓ Touch panel, programmable operation, operating data can be configured according to demands and stored automatically.
- ✓ LED display, easy-to-use interface, simple and convenient operations.
- ✓ Real-time Speed and RCF conversation is suitable for operation.
- ✓ Has a self-locking device, cover safety device, speed safety device and automatic alarm device, etc.
- ✓ Slapstick sealing ring comfort according to GMP. Certified by the US FDA UU



Features:

- ✓ Microcomputer control, brushless drive, programmable operation
- ✓ Metal housing and lid.
- ✓ There is emergency lid release.
- ✓ There are 9 individual acceleration and deceleration stages.
- ✓ Intellectual interface, digital display screen.
- ✓ Stainless steel chamber, anti-causticity force.
- ✓ Equipped with over speed safety device and lid safety device.

Technical Specifications:

Model	YR019-TG
Max Speed	19,000 r/min
Max RCF	24,540g
Max Capacity	750mlx4
Motor	Brushless, Without carbon, DC
Timer Range	0-99min
Noise	<65dB(A)
Power Supply	AC220V 50Hz 5A
Dimension	710×660×490mm (L×W×H)
Net Weight	60kg



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S



KALSTEIN® Centrifuges Tabletop High Speed Centrifuge Model YR019-TG



Optional Accessoires:

Rotor Type	Max Speed (r/min)	Max RCF (xg)	Max Capacity
Angle Rotor	19000	24540	1.5/2.2mlx12
Angle Rotor	16500	19175	5mlx12
Angle Rotor	15000	23120	10mlx12
Angle Rotor	15000	23120	30mlx6
Angle Rotor	13000	17940	50mlx6
Angle Rotor	16500	26050	1.5/2.2mlx24
Angle Rotor	12000	15805	85mlx6
Swing Rotor	10000	11960	100mlx6



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S





 **Kalstein**[®]
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE
ACQUISITION
OF A KALSTEIN
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:

 **Fundación
Jacinto Convit**
<https://www.jacintoconvit.org/>


 **ONE
TREE
PLANTED**

<https://onetreepanted.org>

 **Humatem**
Solidarité et coopération
biomédicale
<https://www.humatem.org/>


**MANIAPURE
FOUNDATION USA**
<https://www.maniapurefoundation.org/>

Articles Kalstein
(Social Networks and Official Websites)



Kalstein Blog
Not just be informed, be fully informed

The 3 articles of the week!

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



1
Read this ARTICLE NOW

How does the climate affect your microscope?

Kalstein Blog

Article of the week:
How does the climate affect your microscope?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



2
Read this ARTICLE NOW

What is the importance of the laboratory refrigerator?

Kalstein Blog

Article of the week:
What is the importance of the laboratory refrigerator?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

Kalstein Training Videos
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Alliances

23 views • 3 years ago



Kalstein's Microscopes

445 views • 1 week ago



PARIS - FRANCE
HEADQUARTER

2 Rue Jean Lantier
75001 Paris - France

Tlf: +33 (0) 01 78 95 87 89

Cel: +33 (0) 6 80 76 07 10

sales@kalstein.eu

<https://kalstein.eu/>
<https://www.kalstein.fr/>

