




LABORATORY LINES

ULTRASONIC CLEANER



Kalstein®

A different accompaniment, at your service





Kalstein®

A different accompaniment, at your service

OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



KALSTEIN® Ultrasonic Cleaner


Ultrasonic Cleaner

Model YR05221



Product Description:

Our desktop ultrasonic cleaner capacity from 2liters to 30liters. Ultrasonic cleaning machine generates ultrasonic frequency electricity, which is converted by energy converter into high-frequency mechanical oscillation that is transferred to cleaning fluid.



Features of Dental Ultrasonic Cleaner:

- ✓ Proprietary MCU-sweep ultrasonic generator driver circuit to present better cleaning performance and to minimize noise level.
- ✓ Proprietary industrial ultrasonic transducer with high Q value and energy conversion efficiency and long service life.
- ✓ Imported high-strength glue, staple-free bonding, and long service life of tank.
- ✓ Digital timing for heating control between 1-99 minutes and temperature setting between 1-80°C.
- ✓ Internal tank fabricated of 1.0mm SUS304 stainless steel to present good performance and long service life.
- ✓ Cleaning basket made of 304 stainless steel and electro-polished.
- ✓ Outer shell made of high-quality stainless steel to present good corrosion resistance.



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
 2 Rue Jean Lantier • 75001 Paris •
 +33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
 KALSTEIN FRANCE, S. A. S



KALSTEIN® Ultrasonic Cleaner Ultrasonic Cleaner



Model	Internal Dimensions	External Dimensions	Volume (L)	Ultrasonic Power(W)	Heating Power(W)	Drain valve	Frequency (KHz)	Time (Min)	Heating Temperature(°C)
YR05220	240*135*100	268*170*290	3.2	120	100	yes			
YR05221	240*135*150	268*170*290	5.0	120	200	yes			
YR05222	300*150*100	325*178*290	4.5	180		yes			
YR05223	300*150*150	325*178*290	6.8	180	200	yes			
YR05224	300*240*150	368*270*290	11	240	300	yes	40	1~99	1~80
YR05225	330*300*150	397*335*317	15	360	400	yes			
YR05226	330*300*200	397*335*368	20	360		yes			
YR05227	500*300*150	560*320*317	22	480	500	yes			
YR05228	500*300*200	560*320*385	30	600	500	yes			



#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S





 **Kalstein**[®]
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE
ACQUISITION
OF A KALSTEIN
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



**Fundación
Jacinto Convit**

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



Humatem

Solidarité et coopération
biomédicale


<https://www.humatem.org/>



**MANIAPURE
FOUNDATION USA**

<https://www.maniapurefoundation.org/>

Articles Kalstein
(Social Networks and Official Websites)



Kalstein Blog
Not just be informed, be fully informed

The 3 articles of the week!

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

1
Read this ARTICLE NOW



How does the climate affect your microscope?

microscope is an instrument that allows observing and measuring objects at a very small scale. This is achieved through a complex system of lenses, which form

Kalstein Blog

Article of the week:
How does the climate affect your microscope?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

2
Read this ARTICLE NOW



What is the importance of the laboratory refrigerator?

refrigerator is a device that is one of the most used equipment in laboratories to maintain, in a controlled environment (refrigerated space), various fluids and substances.

Kalstein Blog

Article of the week:
What is the importance of the laboratory refrigerator?


Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

Kalstein Training Videos
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Alliances

23 views • 3 years ago



Kalstein's Microscopes

445 views • 1 week ago



PARIS - FRANCE
HEADQUARTER

2 Rue Jean Lantier
75001 Paris - France

Tlf: +33 (0) 01 78 95 87 89

Cel: +33 (0) 6 80 76 07 10
sales@kalstein.eu

<https://kalstein.eu/>
<https://www.kalstein.fr/>

