




LABORATORY LINES

WATER DISTILLER



Kalstein®

A different accompaniment, at your service





Kalstein®

A different accompaniment, at your service

OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



#Letsgivemore ♥

Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



Water Distiller Series YR04980



Features:

This device uses the tap-water to produce the pure water and use the electric heating method to produce the distilled water. It is widely used in the schools, laboratories, food processing, biochemistry, agriculture among others.

Characteristics:

- ✓ Main structure of made of stainless steel.
- ✓ Easy operation, steady output
- ✓ Model DZ has special water level control system. It will alarm when the water in evaporation barrel is under the warning water level, and then cut the power automatically, when the water is replenished it continues working.

		YR04981	YR04982	YR04983	YR04984	YR04985	YR04986
Method		water condensation					
Performance	output(L/h)	5	10	20	5	10	20
Structure	evaporator	stainless steel					
	condenser	stainless steel					
	heating device	stainless steel heating tube					
	rated power	4.5kw	7.5kw	15kw	4.5kw	7.5kw	15kw
Controller	automatic alarm	x	x	x	√	√	√
	automatic power cut	x	x	x	√	√	√
	automatic recovery	x	x	x	√	√	√
Safety device	no				water shortage alarm,water shortage power-off,water recovery		
Specification	exterior size (W*D*Hmm)	φ305*730	φ335*830	φ405*980	φ305*730	φ335*830	φ405*980
	packing size (W*D*Hmm)	380*350*810	380*360*900	410*410*1110	380*350*810	380*360*900	410*410*1110
	Power source(50/60Hz)/Rated current	AC220V/20A	AC380V/11A	AC380V/23A	AC220V/20A	AC380V/11A	AC380V/23A
	NW/GW	4.7/6.7kg	6.3/7.8kg	9.5/11kg	5/7kg	7/9kg	10/12kg



#Let'sgive more ♥


Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
Optimum Business Center 450 Rue Baden Powell,
• 34000 Montpellier, France.
Tél: +33 467158849 / +33 680760710/ +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A. S



Articles Kalstein
(Social Networks and Official Websites)



Kalstein Blog
Not just be informed, be fully informed
The 3 articles of the week!
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



1
Read this ARTICLE NOW
How does the climate affect your microscope?
A microscope is an instrument that allows observing and measuring objects that are too small to be seen with the naked eye. This is achieved through a complex system of lenses and mirrors, which form an image of the object being observed.
Kalstein Blog
Article of the week:
How does the climate affect your microscope?
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>




2
Read this ARTICLE NOW
What is the importance of the laboratory refrigerator?
A refrigerator in a laboratory is one of the most important pieces of equipment. It is used to maintain, in a controlled environment (refrigerated space), various fluids and substances.
Kalstein Blog
Article of the week:
What is the importance of the laboratory refrigerator?
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

Kalstein Training Videos
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Alliances

23 views • 3 years ago



Kalstein's Microscopes

445 views • 1 week ago



MONTEPELLIER - FRANCE
SEDE

Optimum Business Center,
450 Rue Baden Powell,
34000 Montpellier - France.

Fax: +33 (0) 4 67 15 88 49

Tel: +33 (0) 6 80 76 07 10

sales@kalstein.eu

<https://kalstein.eu/>
<https://www.kalstein.fr/>

Official Websites for Latin America.

