




MEDICAL LINES

WHEELCHAIRS



Kalstein®

A different accompaniment, at your service





Kalstein®

A different accompaniment, at your service

OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



#Letsgivemore ❤️

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



Electric Wheelchair

Model YR01901X-1



Characteristics:

- ✓ Steel frame
- ✓ Lifting armrest, elevating footrest push handle with rubber handgrip
- ✓ Seat width: 50 cm Seat depth:45 cm Seat height:54 cm Backrest height:75 cm Overall Length:118 cm Overall width: 64cm Overall height:122 cm
- ✓ Rear wheel size: 60cm (same as picture which marked part)
- ✓ Front wheel size: 20cm Folding Height:79 cm Folding Length: 39cm Folding Width:74 cm
- ✓ Gross weight: 50 kg Package Size : 89*41*80cm Lithium battery: 24V 12Ah
- ✓ Load:150KG
- ✓ Warranty: 3 years for frame

- **Relining Backrest**
- **Detachable & Elevating Footrest**



#Letsgivemore 

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S



Manual Wheelchair

Model YR0908LJ



Characteristics:

- ✓ Aluminium frame Flip-up armrest Detachable footrest TPR front wheel
- ✓ Flip-up armrest
- ✓ Detachable footrest
- ✓ TPR Front wheel
- ✓ Handle with brake
- ✓ Solid castor, PU rear wheel handle with brake
- ✓ PUSH HANDLE WITH RUBBER HANDGRIP
- ✓ Seat width: 50 cm
- ✓ Seat depth: 41 cm
- ✓ Seat height: 48 cm
- ✓ Backrest height: 38cm
- ✓ Overall Length: 104 cm
- ✓ Front wheel: 20cm
- ✓ Rear wheel: 60cm
- ✓ Load: 150kg
- ✓ Packing size: 88 * 27 * 82cm Warranty: 3 years for frame



12KG/14.8KG	100KG	104cm	61cm	88cm	46cm	41cm	6 inches	22 inches	38cm	48cm	86*80*26cm	82*27*88cm



#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
 2 Rue Jean Lantier • 75001 Paris •
 +33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
 KALSTEIN FRANCE, S. A. S



Kids Wheelchair

Model YR0M2207



Characteristics:

- ✓ Steel frame Flip-up armrest
- ✓ Detachable footrest
- ✓ PU rubber front wheel
- ✓ PU rubber rear wheel handle with brake
- ✓ PUSH HANDLE WITH RUBBER HANDGRIP
- ✓ Seat width: 38 cm
- ✓ Seat depth: 45 cm
- ✓ Seat height: 45 cm
- ✓ Backrest height: 45cm
- ✓ Overall Length: 104 cm
- ✓ Front wheel: 20cm Rear wheel: 51cm
- ✓ Load: 100kg
- ✓ Packing size: 82*24*87 cm Warranty: 3 years for frame



- **Flip-up armrest**
- **Multicolours for frame**
- **Detachable footrest**



#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S





 **Kalstein**[®]
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE
ACQUISITION
OF A KALSTEIN
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



**Fundación
Jacinto Convit**

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



Humatem

Solidarité et coopération
biomédicale


<https://www.humatem.org/>



**MANIAPURE
FOUNDATION USA**

<https://www.maniapurefoundation.org/>

Articles Kalstein
(Social Networks and Official Websites)



Kalstein Blog
Not just be informed, be fully informed
The 3 articles of the week!
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



1
Read this ARTICLE NOW
How does the climate affect your microscope?
A microscope is an instrument that allows observing and measuring objects that are too small to be seen with the naked eye. This is achieved through a complex system of lenses and mirrors, which form...

Kalstein Blog
Article of the week:
How does the climate affect your microscope?
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



2
Read this ARTICLE NOW
What is the importance of the laboratory refrigerator?
A refrigerator is a device that is one of the most common pieces of equipment used in laboratories. It is a controlled environment (refrigerated space), various fluids and substances...

Kalstein Blog
Article of the week:
What is the importance of the laboratory refrigerator?
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

Kalstein Training Videos
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Alliances
23 views • 3 years ago



Kalstein's Microscopes
445 views • 1 week ago



PARIS - FRANCE
HEADQUARTER
2 Rue Jean Lantier
75001 Paris - France

Tlf: +33 (0) 01 78 95 87 89
Cel: +33 (0) 6 80 76 07 10
sales@kalstein.eu

<https://kalstein.eu/>
<https://www.kalstein.fr/>

